Division of Environment Kansas Health and Environmental Laboratories 6810 SE Dwight Street Topeka, KS 66620



Phone: 785-296-1620 Fax: 785-559-5205 www.kdheks.gov/labs

## **REQUISITION FOR LABORATORY SPECIMEN KITS and SUPPLIES**

- Refer to the Manual of Laboratory Tests or call (785) 296-1620 regarding specimen submission
- Return Completed form by fax or email to: (785) 559-5205 or KDHE.KHEL Help@ks.gov

Enter the quantity needed on the line next to the item

Health Specimen Submission Forms	
Universal Specimen Submission Forms	Neonatal Brochure
Neonatal Screening	English
Phenylketonuria Monitoring	Spanish
Health	Specimen Kits
Virology/Serology	Tuberculosis
Multi-tube Container with Mailing Box (Blood – Ambient Valid for Hepatitis ONLY)	TB Mailer with Sputum Collection tube QuantiFERON (QFT tubes only) Multi-tube Container with Mailing Box (Blood – Ambient, Valid for QFT)
Multi-tube Container with Mailing Box (Chlamydia/Gonorrhea – Ambient)	
Mailing Cooler (Cold shipper)  (Serum – Cold, Valid for ALL Serology)	Bacterial Enteric Mailer with Cary Blair (Category B)
Monkeypox Collection Kit (Swabs and empty tubes)	Bacterial Isolate (Category B) Bacterial Isolate (Category A / IDS)
Other H	lealth Supplies
EDTA Purple Top Blood Tubes (Venous Blood Metals)  K2 EDTA Microtainer (Capillary Blood Metals)  Serum Collection SST tubes (Gold-top)  Screw-cap Serum Tubes (Pour-off tube)  Viral Transport Medium (VTM)  Other (Please specify)	Miscellaneous Category B Shipper ambient (box only)  Miscellaneous Category B Shipper cold (box & cooler only)  Lancet (Capillary Blood Lead collection only – not NBS)  Specimen Transport Bags  Nasopharyngeal Swab
Facility ID Number: Facility Name:	E-mail completed form to KDHE.KHEL_Help@ks.gov
Attention: Address: State: Zip: Telephone:	Or Fax to 785-559-5205